

Anti-ZMPSTE24 antibody (450-480) {Biotin} (STJ503599)
STJ503599**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZMPSTE24 (450-480) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	230709
Gene Symbol	Zmpste24
Uniprot ID	FACE1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 450-480 on mouse ZMPSTE24 protein.
Immunogen Region	450-480
Specificity	
Immunogen Sequence	