

5HT1A Blocking Peptide for STJ500001 peptide (STJ503611)

STJ503611

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	5HT1A Blocking Peptide for STJ500001 is synthetically produced from the 200-250 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	24473
Gene Symbol	Htr1a
Uniprot ID	5HT1A_RAT
Immunogen	Synthetic peptide corresponding to amino acids within region 200-250 on rat 5-hydroxytryptamine receptor 1A protein.
Immunogen Region	200-250
Specificity	This blocking peptide is recommended for use in combination with 5HT1A antibody, STJ500001
Immunogen Sequence	