

Alpha 1 Fetoprotein Blocking Peptide peptide (STJ503631)

STJ503631

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Alpha 1 Fetoprotein Blocking Peptide is synthetically produced from the Full length sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	174
Gene Symbol	AFP
Uniprot ID	FETA_HUMAN
Immunogen	Full length native protein (purified) corresponding to Human alpha 1 Fetoprotein. Purified human alpha 1 Fetoprotein.
Immunogen Region	Full length
Specificity	
Immunogen Sequence	