

Adenine Nucleotide Translocator 2 Positive Control peptide (STJ503647)

STJ503647

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Adenine Nucleotide Translocator 2 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	292
Gene Symbol	SLC25A5
Uniprot ID	ADT2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	