

ARF1 Blocking Peptide for STJ500153 peptide (STJ503683)

STJ503683

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ARF1 Blocking Peptide for STJ500153 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64310
Gene Symbol	Arf1
Uniprot ID	ARF1_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Rat ARF1 protein. The peptide was covalently modified post-synthetically to achieve desired antigenicity.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with ARF1 antibody, STJ500153
Immunogen Sequence	