

ARHGEF7 Positive Control for STJ500159 peptide (STJ503702)

STJ503702

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ARHGEF7 Positive Control for STJ500159 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8874
Gene Symbol	ARHGEF7
Uniprot ID	ARHG7_HUMAN
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with ARHGEF7 antibody STJ500159.
Immunogen Sequence	