

Bestrophin 3 Blocking Peptide peptide (STJ503721)

STJ503721

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Bestrophin 3 Blocking Peptide is synthetically produced from the C-Term sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	382427
Gene Symbol	Best3
Uniprot ID	BEST3_MOUSE
Immunogen	Synthetic peptide taken within the last 50 amino acids from the C-terminus of rat Bestrophin-3 protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	