

BMP11 Positive Control peptide (STJ503749)

STJ503749

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	BMP11 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10220
Gene Symbol	GDF11
Uniprot ID	GDF11_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	