

CD36 Blocking Peptide for STJ500464 peptide (STJ503833)

STJ503833

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	CD36 Blocking Peptide for STJ500464 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	948
Gene Symbol	CD36
Uniprot ID	CD36_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human CD36 protein.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with CD36 antibody, STJ500464
Immunogen Sequence	