

CXCR1 Blocking Peptide for STJ500645 peptide (STJ503879)

STJ503879

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	CXCR1 Blocking Peptide for STJ500645 is synthetically produced from the 272-350 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3577
Gene Symbol	CXCR1
Uniprot ID	CXCR1_HUMAN
Immunogen	Synthetic peptide, corresponding to a region within C terminal amino acids 272-350 of Human CXCR1.
Immunogen Region	272-350
Specificity	This blocking peptide is recommended for use in combination with CXCR1 antibody, STJ500645
Immunogen Sequence	