

Dopamine D2 Receptor Blocking Peptide peptide (STJ503942)

STJ503942

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Dopamine D2 Receptor Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1813
Gene Symbol	DRD2
Uniprot ID	DRD2_HUMAN
Immunogen	Synthetic peptide corresponding to unique sequence on human Dopamine D2 Receptor protein.
Immunogen Region	
Specificity	
Immunogen Sequence	