

Dopamine Receptor D4 Blocking Peptide for STJ500717 peptide (STJ503944)

STJ503944

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Dopamine Receptor D4 Blocking Peptide for STJ500717 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	13491
Gene Symbol	Drd4
Uniprot ID	DRD4_MOUSE
Immunogen	Synthetic peptide corresponding to unique sequence on Dopamine Receptor D4 protein.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with Dopamine Receptor D4 antibody, STJ500717
Immunogen Sequence	