

DCC Blocking Peptide for STJ500740 peptide (STJ503948)

STJ503948

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	DCC Blocking Peptide for STJ500740 is synthetically produced from the 1330-1380 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1630
Gene Symbol	DCC
Uniprot ID	DCC_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1330-1380 on rat DCC protein.
Immunogen Region	1330-1380
Specificity	This blocking peptide is recommended for use in combination with DCC antibody, STJ500740
Immunogen Sequence	