

DOC2B Blocking Peptide for STJ500786 peptide (STJ503960)

STJ503960

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	DOC2B Blocking Peptide for STJ500786 is synthetically produced from the 360-412 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	81820
Gene Symbol	Doc2b
Uniprot ID	DOC2B_RAT
Immunogen	Synthetic Peptide taken within amino acid region 360-412 on DOC2b gene.
Immunogen Region	360-412
Specificity	This blocking peptide is recommended for use in combination with DOC2B antibody, STJ500786
Immunogen Sequence	