

Dopamine Receptor D1 Positive Control for STJ500706 peptide (STJ503964)

STJ503964

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Dopamine Receptor D1 Positive Control for STJ500706 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	Drd1
Uniprot ID	DRD1_RAT
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with Dopamine Receptor D1 antibody STJ500706.
Immunogen Sequence	