

DOC2B Positive Control for STJ500788 peptide (STJ503980)

STJ503980

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	DOC2B Positive Control for STJ500788 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	81820
Gene Symbol	Doc2b
Uniprot ID	DOC2B_RAT
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with DOC2B antibody STJ500788.
Immunogen Sequence	