

E-Cadherin Blocking Peptide for STJ500839 peptide (STJ503992)

STJ503992

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	E-Cadherin Blocking Peptide for STJ500839 is synthetically produced from the 50-100 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human E-cadherin 1 protein.
Immunogen Region	50-100
Specificity	This blocking peptide is recommended for use in combination with E-Cadherin antibody, STJ500839
Immunogen Sequence	