

ELOVL4 Blocking Peptide for STJ500878 peptide (STJ504004)

STJ504004

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ELOVL4 Blocking Peptide for STJ500878 is synthetically produced from the 294-301 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6785
Gene Symbol	ELOVL4
Uniprot ID	ELOV4_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 294-301 of Human ELOVL4. Sequence: SEKQLMIENGKKQKNGKAKGDC
Immunogen Region	294-301
Specificity	This blocking peptide is recommended for use in combination with ELOVL4 antibody, STJ500878
Immunogen Sequence	