

SLC29A3 Blocking Peptide for STJ500885 peptide (STJ504005)

STJ504005

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	SLC29A3 Blocking Peptide for STJ500885 is synthetically produced from the 1-60 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55315
Gene Symbol	SLC29A3
Uniprot ID	S29A3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-60 on ENT3 protein.
Immunogen Region	1-60
Specificity	This blocking peptide is recommended for use in combination with SLC29A3 antibody, STJ500885
Immunogen Sequence	