

ERG Blocking Peptide for STJ500948 peptide (STJ504021)

STJ504021

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ERG Blocking Peptide for STJ500948 is synthetically produced from the 400-450 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2078
Gene Symbol	ERG
Uniprot ID	ERG_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Transcriptional regulator ERG protein.
Immunogen Region	400-450
Specificity	This blocking peptide is recommended for use in combination with ERG antibody, STJ500948
Immunogen Sequence	