

ELOVL4 Positive Control for STJ500874 peptide (STJ504041)

STJ504041

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Biomolecules |
| Short Description | ELOVL4 Positive Control for STJ500874 is synthetically produced from the sequence and is suitable for use in western blot applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | WB: 1:500 |
| Formulation | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 6785 |
| Gene Symbol | ELOVL4 |
| Uniprot ID | ELOV4_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | This is positive control is recommended for use in combination with ELOVL4 antibody STJ500874. |
| Immunogen Sequence | |