

FAM38A Blocking Peptide peptide (STJ504063)

STJ504063

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	FAM38A Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9780
Gene Symbol	PIEZO1
Uniprot ID	PIEZ1_HUMAN
Immunogen	Synthetic peptide taken between TMD14 and TMD15 on human FAM38A protein.
Immunogen Region	
Specificity	
Immunogen Sequence	