

FEL D1 Major allergen 1 polypeptide chain 1 Blocking Peptide for STJ500997 peptide (STJ504067)

STJ504067

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	FEL D1 Major allergen 1 polypeptide chain 1 Blocking Peptide for STJ500997 is synthetically produced from the 1-50 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1-50 on feline Major cat allergen I polypeptide chain 2.
Immunogen Region	1-50
Specificity	This blocking peptide is recommended for use in combination with FEL D1 Major allergen 1 polypeptide chain 1 antibody, STJ500997
Immunogen Sequence	