

GPCR GPR120 FFAR4-412AP Blocking Peptide peptide (STJ504074)

STJ504074

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	GPCR GPR120 FFAR4-412AP Blocking Peptide is synthetically produced from the 220-300 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	338557
Gene Symbol	FFAR4
Uniprot ID	FFAR4_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 220-300 on human FFAR4 protein.
Immunogen Region	220-300
Specificity	
Immunogen Sequence	