

C11orf1 FLJ-112AP Blocking Peptide peptide (STJ504083)

STJ504083

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	C11orf1 FLJ-112AP Blocking Peptide is synthetically produced from the 50-100 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64776
Gene Symbol	CFAP68
Uniprot ID	CFA68_HUMAN
Immunogen	Mid-region epitope synthetic peptide taken within amino acid region 50-100 on human C11orf1 protein
Immunogen Region	50-100
Specificity	
Immunogen Sequence	