

FOXP3 Blocking Peptide peptide (STJ504085)

STJ504085

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	FOXP3 Blocking Peptide is synthetically produced from the 229-250 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64776
Gene Symbol	CFAP68
Uniprot ID	CFA68_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 229-250 of Mouse Forkhead Box Protein 3 (winged helix transcription factor 3 protein). The peptide was post synthetically modified to achieve desired antigenicity before coupling to a carrier protein. Se
Immunogen	229-250
Region	
Specificity	
Immunogen	
Sequence	