

FUBP1 Blocking Peptide for STJ501094 peptide (STJ504088)

STJ504088

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	FUBP1 Blocking Peptide for STJ501094 is synthetically produced from the 600-644 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8880
Gene Symbol	FUBP1
Uniprot ID	FUBP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 600-644 on human FUBP1 protein.
Immunogen Region	600-644
Specificity	This blocking peptide is recommended for use in combination with FUBP1 antibody, STJ501094
Immunogen Sequence	