

DDIT3 Blocking Peptide for STJ501131 peptide (STJ504122)

STJ504122

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	DDIT3 Blocking Peptide for STJ501131 is synthetically produced from the 60-160 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1649
Gene Symbol	DDIT3
Uniprot ID	DDIT3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 60-160 on human GADD153 protein
Immunogen Region	60-160
Specificity	This blocking peptide is recommended for use in combination with DDIT3 antibody, STJ501131
Immunogen Sequence	