

GLUT12 Blocking Peptide for STJ501202 peptide (STJ504140)

STJ504140

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	GLUT12 Blocking Peptide for STJ501202 is synthetically produced from the C-Term sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide from the C-terminus of Rat GLUT12.
Immunogen	C-Term
Region	
Specificity	This blocking peptide is recommended for use in combination with GLUT12 antibody, STJ501202
Immunogen	
Sequence	