

## Golgin A1 Blocking Peptide for STJ501232 peptide (STJ504150)

STJ504150

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Golgin A1 Blocking Peptide for STJ501232 is synthetically produced from the 650-700 sequence and is suitable for use in western blot applications.
<b>Applications</b>	Immunodepletion/Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">76899</a>
<b>Gene Symbol</b>	<a href="#">Golga1</a>
<b>Uniprot ID</b>	<a href="#">GOGA1_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 650-700 on mouse Golgin subfamily A member 1 isoform b protein.
<b>Immunogen Region</b>	650-700
<b>Specificity</b>	This blocking peptide is recommended for use in combination with Golgin A1 antibody, STJ501232
<b>Immunogen Sequence</b>	