

Golgin A3 Blocking Peptide peptide (STJ504152)

STJ504152

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Golgin A3 Blocking Peptide is synthetically produced from the 1400-1450 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2802
Gene Symbol	GOLGA3
Uniprot ID	GOGA3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1400-1450 on human Golgin subfamily A member 3 isoform 1 protein.
Immunogen Region	1400-1450
Specificity	
Immunogen Sequence	