

Glucose Transporter GLUT2 Positive Control peptide (STJ504193)

STJ504193

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Glucose Transporter GLUT2 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25139
Gene Symbol	Slc2a4
Uniprot ID	GLUT4_RAT
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	