

HAND2 Blocking Peptide peptide (STJ504216)

STJ504216

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	HAND2 Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64637
Gene Symbol	Hand2
Uniprot ID	HAND2_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on Heart and neuronal crest derivative protein 2.
Immunogen Region	
Specificity	
Immunogen Sequence	