

Hepatitis C virus Core NS3 + NS4 Positive Control for STJ501321 peptide (STJ504235)

STJ504235

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Hepatitis C virus Core NS3 + NS4 Positive Control for STJ501321 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with Hepatitis C virus Core NS3 + NS4 antibody STJ501321.
Immunogen Sequence	