

Outer membrane protein 26 H pylori Blocking Peptide for STJ501362 peptide (STJ504237)

STJ504237

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Outer membrane protein 26 H pylori Blocking Peptide for STJ501362 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Heat killed and denatured Helicobacter pylori extracts
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with Outer membrane protein 26 H. pylori antibody, STJ501362
Immunogen Sequence	