

IL17RB Blocking Peptide for STJ501428 peptide (STJ504251)

STJ504251

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	IL17RB Blocking Peptide for STJ501428 is synthetically produced from the 300-400 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55540
Gene Symbol	IL17RB
Uniprot ID	I17RB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-400 on IL17RB protein.
Immunogen Region	300-400
Specificity	This blocking peptide is recommended for use in combination with IL17RB antibody, STJ501428
Immunogen Sequence	