

ITK Blocking Peptide for STJ501511 peptide (STJ504271)

STJ504271

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ITK Blocking Peptide for STJ501511 is synthetically produced from the 550-600 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3702
Gene Symbol	ITK
Uniprot ID	ITK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 550-600 on human ITK protein.
Immunogen Region	550-600
Specificity	This blocking peptide is recommended for use in combination with ITK antibody, STJ501511
Immunogen Sequence	