

KLF6 isoform C Blocking Peptide for STJ501558 peptide (STJ504286)

STJ504286

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Biomolecules |
| Short Description | KLF6 isoform C Blocking Peptide for STJ501558 is synthetically produced from the 200-250 sequence and is suitable for use in western blot applications. |
| Applications | Immunodepletion/Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Synthetic peptide taken within amino acid region 200-250 on human KLF6 protein. |
| Immunogen Region | 200-250 |
| Specificity | This blocking peptide is recommended for use in combination with KLF6 isoform C antibody, STJ501558 |
| Immunogen Sequence | |