

MAP3K10 Blocking Peptide for STJ501665 peptide (STJ504358)

STJ504358

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MAP3K10 Blocking Peptide for STJ501665 is synthetically produced from the 900-950 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4294
Gene Symbol	MAP3K10
Uniprot ID	M3K10_HUMAN
Immunogen	Synthetic peptide corresponding to the C-terminal taken within amino acid region 900-950 on human MAP3K10 protein.
Immunogen Region	900-950
Specificity	This blocking peptide is recommended for use in combination with MAP3K10 antibody, STJ501665
Immunogen Sequence	