

MAP3K12 Blocking Peptide for STJ501666 peptide (STJ504359)

STJ504359

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MAP3K12 Blocking Peptide for STJ501666 is synthetically produced from the 270-320 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7786
Gene Symbol	MAP3K12
Uniprot ID	M3K12_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 270-320 on MAP3K12 protein
Immunogen Region	270-320
Specificity	This blocking peptide is recommended for use in combination with MAP3K12 antibody, STJ501666
Immunogen Sequence	