

MC3 Receptor Blocking Peptide peptide (STJ504363)

STJ504363

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MC3 Receptor Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4157
Gene Symbol	MC1R
Uniprot ID	MSHR_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on MC3r protein.
Immunogen Region	
Specificity	
Immunogen Sequence	