

MRAP2 Blocking Peptide for STJ501796 peptide (STJ504383)

STJ504383

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MRAP2 Blocking Peptide for STJ501796 is synthetically produced from the 70-120 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	77037
Gene Symbol	Mrap
Uniprot ID	MRAP_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 70-120 on mouse Melanocortin-2 receptor accessory protein.
Immunogen Region	70-120
Specificity	This blocking peptide is recommended for use in combination with MRAP2 antibody, STJ501796
Immunogen Sequence	