

MRP1 Blocking Peptide for STJ501801 peptide (STJ504385)

STJ504385

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	MRP1 Blocking Peptide for STJ501801 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	24565
Gene Symbol	Abcc1
Uniprot ID	MRP1_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on rat multidrug resistance-associated protein 1.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with MRP1 antibody, STJ501801
Immunogen Sequence	