

Myocilin Blocking Peptide peptide (STJ504396)

STJ504396

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Myocilin Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4653
Gene Symbol	MYOC
Uniprot ID	MYOC_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Myocilin protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	