

## NME2 + NM23A Blocking Peptide peptide (STJ504404)

STJ504404

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	NME2 + NM23A Blocking Peptide is synthetically produced from the 100-150 sequence and is suitable for use in western blot applications.
<b>Applications</b>	Immunodepletion/Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">4830</a>
<b>Gene Symbol</b>	<a href="#">NME1</a>
<b>Uniprot ID</b>	<a href="#">NDKA_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 100-150 on NM23A protein selective for both NDKA and NDKB.
<b>Immunogen Region</b>	100-150
<b>Specificity</b>	
<b>Immunogen Sequence</b>	