

**MRGX1 Positive Control for STJ501799 peptide (STJ504459)**

STJ504459

**GENERAL INFORMATION**

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	MRGX1 Positive Control for STJ501799 is synthetically produced from the sequence and is suitable for use in western blot applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: 1:500
<b>Formulation</b>	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	282547
<b>Gene Symbol</b>	Mrgprx1
<b>Uniprot ID</b>	MRGX1_RAT
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	This is positive control is recommended for use in combination with MRGX1 antibody STJ501799.
<b>Immunogen Sequence</b>	