

ADCY4 Blocking Peptide for STJ502069 peptide (STJ504515)

STJ504515

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ADCY4 Blocking Peptide for STJ502069 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on rat Particulate Adenylyl Cyclase 4 protein.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with ADCY4 antibody, STJ502069
Immunogen Sequence	