

ADCY6 Blocking Peptide for STJ502072 peptide (STJ504517)

STJ504517

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ADCY6 Blocking Peptide for STJ502072 is synthetically produced from the 13-27 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25289
Gene Symbol	Adcy6
Uniprot ID	ADCY6_RAT
Immunogen	Synthetic peptide corresponding to N-terminal amino acids 13-27 of rat Particulate Adenylyl Cyclase 6 protein. Peptide was covalently modified and conjugated to KLH. Sequence: DERKTAWGERNGQKR
Immunogen Region	13-27
Specificity	This blocking peptide is recommended for use in combination with ADCY6 antibody, STJ502072
Immunogen Sequence	