

AKT1 Blocking Peptide peptide (Phospho) (STJ504523)

STJ504523

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-AKT1 Blocking Peptide is synthetically produced from the 430-480 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	207
Gene Symbol	AKT1
Uniprot ID	AKT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 430-480 on human Rac protein kinase-alpha protein.
Immunogen Region	430-480
Specificity	
Immunogen	
Sequence	