

Annexin A2 Blocking Peptide peptide (Phospho) (STJ504527)

STJ504527

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Annexin A2 Blocking Peptide is synthetically produced from the 1-50 aa sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	302
Gene Symbol	ANXA2
Uniprot ID	ANXA2_HUMAN
Immunogen	Synthetic peptide corresponding to Human Annexin A2 protein with phosphorylated serine 12 and threonine 18 residues. Synthetic peptide was taken within amino acid region 1-50.
Immunogen Region	1-50 aa
Specificity	
Immunogen Sequence	